

IHE Integration Statement (1/3)		Date	23 Feb. 2015
Vendor	Product Name	Version	
PSP Corporation.	EV Insite net	3.1.1.200	
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow (SWF)	Image Manager/Image Archive (IM/IA)	(Japanese Character Set)	
	Performed Procedure Step Manager (PPSM)	No options defined	
	Image Display (ID)	No options defined	
Scheduled Workflow .b(SWF.b)	Image Manager/Image Archive (IM/IA)	No options defined	
Patient Information Reconciliation (PIR)	Image Manager/Image Archive (IM/IA)	No options defined	
	Performed Procedure Step Manager (PPSM)	No options defined	
Consistent Presentation of Images (CPI)	Image Manager/Image Archive (IM/IA)	No options defined	
	Image Display (ID)	No options defined	
	Evidence Creator(EC)	No options defined	
Access to Radiology Information (ARI)	Image Manager/Image Archive (IM/IA)	No options defined	
	Image Display (ID)	No options defined	
Key Image Note(KIN)	Image Manager/Image Archive (IM/IA)	No options defined	
	Image Display (ID)	No options defined	
	Evidence Creator(EC)	No options defined	
NM Image (NM)	Image Manager/Image Archive (IM/IA)	No options defined	

IHE Integration Statement (2/3)		Date	23 Feb. 2015
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Mammography Image (MAMMO)	Image Manager/Image Archive (IM/IA)	No options defined	
	Image Display (ID)	none	
	Print Composer (PC)	No options defined	
Portable Data for Imaging (PDI)	Portable Media Creator (PMC)	Web Content	
	Portable Media Importer (PMI)	No options defined	
	Image Display (ID)	No options defined	
	Print Composer (PC)	No options defined	
	Display (DISPLAY)	No options defined	
Import Reconciliation Workflow Profile (IRWF)	Image Manager/Image Archive (IM/IA)	No options defined	
	Importer (IMPORTER)	No options defined	
	Performed Procedure Step Manager (PPSM)	No options defined	
Radiation Exposure Monitoring (REM)	Image Manager/Image Archive (IM/IA)	No options defined	
Basic Radiation Therapy Objects (BRTO)	Archive	None	
Multimodality Image Registration for Radiation Oncology (MMRO)	Archive	None	
Audit Trail and Node Authentication (ATNA)	Audit Recode Repository	None	
	Secure Application	None	

IHE Integration Statement (3/3)		Date	23 Feb. 2015
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Cardiac Catheterization Workflow (CATH)	Image Manager/Image Archive (IM/IA)	None	
	Image Display (ID)	None	
	Performed Procedure Step Manager (PPSM)	None	
Echocardiography Workflow (ECHO)	Image Manager/Image Archive (IM/IA)	Intermittently Connected Modality	
	Image Display (ID)	None	
	Performed Procedure Step Manager (PPSM)	None	
Intravascular Imaging Workflow (IVI)	Image Manager/Image Archive (IM/IA)	Intermittently Connected Modality	
	Performed Procedure Step Manager (PPSM)	No options defined	
	Image Display (ID)	None	
ED-CARD	Image Manager/Image Archive (IM/IA)	None	
	Performed Procedure Step Manager (PPSM)	None	
	Image Display (ID)	None	
STRESS	Image Manager/Image Archive (IM/IA)	None	
	Performed Procedure Step Manager (PPSM)	None	
	Image Display (ID)	None	
ECG	Display(DISPLAY)	None	
Consistent Time (CT)	Time Client (TC)	None	
Internet address for vendor's IHE information: http://www.psp.co.jp			
Links to Standards Conformance Statements for the Implementation			
HL7	None		
DICOM	http://www.psp.co.jp/download.html		
Links to general information on IHE			
In North America: http://www.ihe.net		In Europe: http://www.ihe-europe.org	In Japan: http://www.ihe-j.org
Remarks:			