

IHE Integration Statement (1/1)		Date	23 Feb. 2015
Vendor	Product Name	Version	
PSP Corporation.	ARISStation	2.50.395.250	
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Import Reconciliation Workflow Profile (IRWF)	Order Filler(OF)	none	
	Performed Procedure Step Manager (PPSM)	No options defined	
Cardiac Catheterization Workflow (CATH)	Order Filler (OF)	none	
	Performed Procedure Step Manager (PPSM)	No options defined	
Echocardiography Workflow (ECHO)	Order Filler (OF)	none	
	Performed Procedure Step Manager (PPSM)	No options defined	
STRESS	Order Filler (OF)	none	
	Performed Procedure Step Manager (PPSM)	No options defined	
Scheduled Workflow (SWF)	Performed Procedure Step Manager (PPSM)	No options defined	
	Order Filler (OF)	none	
Scheduled Workflow .b(SWF.b)	Order Filler (OF)	none	
Consistent Time (CT)	Time Client (TC)	none	
Internet address for vendor's IHE information: http://www.psp.co.jp			
Links to Standards Conformance Statements for the Implementation			
HL7	None		
DICOM	http://www.psp.co.jp/download.html		
Links to general information on IHE			
In North America: http://www.ihe.net		In Europe: http://www.ihe-europe.org	In Japan: http://www.ihe-j.org
Remarks:			